

6-Drawer Tall Anesthesia Cart (UTGSU-333369-DKB)

 waterloohealthcare.com/product/6-drawer-tall-utgsu-333369-dkb



Waterloo anesthesia carts have been widely used in the OR for decades. With a variety of organizing accessories and drawer trays, anesthesiologists have critical supplies at their fingertips.

KEY FEATURES:

- Variety of locking options available including dual-credential/auto-relocking controlled substance drawer
- Single pull out shelf
- Economical steel construction


SKU: UTGSU-333369-DKB Category: Anesthesia Carts: The Best Carts for Your Medical Facility

SPECIFICATIONS

Frame Type: STEEL
 Frame Size: TALL
 Dimensions: 29"w x 24.5"d x 42"h



Drawer Dimensions: 22"w x 16.5"d



Drawer Height: 3-3-3-3-6-9

3" 4 drawers

6" 1 drawer

9" 1 drawer

12"


Frame Color: LIGHT GRAY




Drawer Color: DARK BLUE



Lock System: PUSH BUTTON




Top: REMOVABLE PLASTIC



Bumpers: STABILIZER BASE



Casters: 5" STANDARD



Compare Security Locking Systems

GETTING THE RIGHT ACCESS METHOD	Key Lock	Lockbox	Mechanical Push Button	Electronic Keypad	WiFi Enabled Keypad	WiFi Enabled Keypad with Controlled Substance Drawer
	\$	\$	\$	\$	\$	\$
Standard Key Access	●					
Tamper Evidence		●				
Individual Drawer Locking with Key or Breakaway Seal	●	●	●	●	●	●
Keyless		●	●	●	●	●
Automatic Relocking				●	●	●
Proximity/Employee Badge Access				●	●	●
Ability to Track Lock/Unlock Events				●	●	●
Ability to Control User Access to Specific Times During the Day				●	●	●
Updates to Lock Performed MANUALLY for use of keypad				●		
Updates to Lock Performed AUTOMATICALLY via WiFi					●	●
Maximum Unique User Count				250	3000	3000
Requires LockManager Software				●	●	●

● Standard Feature ● Optional Feature

[Click Here to Learn More About Our Locking Systems >>](#)